



MediZysis India Pvt. Ltd.

one stop solution

Steriguard I.C

QAC BASED DISINFECTANT AND DEODORISER, CLEANER FOR AIR AND SURFACE

Composition

Benzalkonium chloride IP - 0.5% w/v
Isopropyl Alcohol IP(2-propenol) IP - 70% w/v
Perfume

Steriguard I.C

Steriguard IC is a quaternary ammonium & alcohol based instant disinfectant suitable for electronic devices and all other equipments where water based products can not be used. It is Broad spectrum, Bactericidal, Fungicidal, Virucidal, cleaner



Features

- ◆ Aldehyde free, Eco-friendly
- ◆ Easy to use
- ◆ Long Residual Effect
- ◆ Potent, synergistic formulation
- ◆ Fast action
- ◆ Non- corrosive



MediZysis India Pvt. Ltd.
transforming hygiene...

Customer Care: 011-22343535, 9425628989

Email: medizysis@gmail.com

Application

STERIGUARD IC is highly effective against wide range of microorganisms and ready to use liquid disinfectant. It cleans and removes dust, mold, other common soils very easily and quickly from hard and non porous surfaces

STERIGUARD IC is ideal for OT equipment's, Dental Chairs, Hospital beds, Monitors, NICU pumps, Stethoscopes, Incubators, OT table, OT lamp, and pharmaceutical sterile sections.

Direction of Use

- ◆ **Surface** - Spray **STERIGUARD I.C** on the surface to be disinfected and wipe with sterile cloth. It has got prolonged residual effect which insures good results against all the microorganisms.
- ◆ **Environment** - Spray **STERIGUARD I.C** to contaminated rooms, clinical chambers, pathological laboratories with the help of a sprayer

Contact time

0-60 seconds for invisible soiled areas.
15 minutes for visible soiled areas

Precautions

For External Use only.
Avoid direct contact with eyes & skin.
Keep away from fire or flame.
Seek medical help in case of emergency.
Keep out of reach of children.
Highly Flammable.

Physical & Chemical Properties

Colour	-	Colourless
PH	-	8-9
Odour	-	Alcoholic
Storage	-	At room Temperature

Material Compatibility

STERIGUARD I.C is compatible with most of the metals, Plastics elastomers and adhesives.



MediZysis India Pvt. Ltd.
transforming hyziene...

B-1/209, Harsh Vihar, Delhi-110093
Customer Care: 011-22343535, 9425628989
Email: medizysis@gmail.com

Name of Product - STERIGUARD I.C
Pack Size - 500 ML